

MorphoCatcher: a multiple-alignment based web tool for target selection and designing taxon-specific primers in the loop-mediated isothermal amplification method

Fedor V. Shirshikov^{1,*}, Yuri A. Pekov² and Konstantin A. Miroshnikov¹

¹Shemyakin & Ovchinnikov Institute of Bioorganic Chemistry RAS, Moscow, Russia

²Lomonosov Moscow State University, Moscow, Russia

*To whom correspondence should be addressed.

Contact: shrshkv@ya.ru

Introduction

In this short manual we review two case studies of using the MorphoCatcher service for development of taxon-specific or group-specific primer sets and detection of the following potato pathogens by the loop-mediated isothermal amplification (LAMP) method:

- Bacterial potato pathogens from the *Dickeya* genus using an orthologous gene;
- Potato virus Y (PVY) using a whole-genome sequence alignment.

The MorphoCatcher protocol deals with well-established standards of the sequence file formatting. The user can handle FASTA-files with small modifications in the header. The same modifications are needed during formatting of multiple sequence alignment (MSA) using the Clustal Omega web service, because of most deposited nucleotide sequences contain the FASTA-header with spaces that cause trimming of the header after alignment.

It is necessary to understand that the case studies mostly need to illustrate main opportunities of the software developed. We agree that there are more suitable targets for diagnostics of the target taxon. Here, the correct formatting of text files for data processing also shown for each step of the protocol.


There is no efficient bioinformatics tool for prediction of good primer set. Therefore, we recommend to test different primer sets in laboratory, because at present this is a best way to validate specificity and sensitivity of your LAMP assay. In some cases, the column purification of inner primers is needed to enhance sensitivity, rapidity, and reproducibility of the amplification.

We hope that the MorphoCatcher service will help to our colleagues to create new diagnostic assays not only for the LAMP method, but also for other perspective amplification methods. Please, feel free to contact us by e-mail with any questions and feedback. Together we can make this web tool better for the molecular diagnostics community.

A command line version of the MorphoCatcher tool and set of tutorial files are available at the GitHub public repository: <https://github.com/shrshkv/MorphoCatcher>.

Case study 1: the bacterial pathogens and orthologous target gene

To perform the primer design correctly, we recommend to follow these steps:



Step 1. Defining closely related bacterial species.

One of the most widespread and dangerous bacterial pathogens of potato belong to the *Pectobacteriaceae* family [1]. The same potato tubers can be contaminated by several species simultaneously, for example by *Dickeya* and *Pectobacterium* spp. Therefore, we should to include some species from the abovementioned genera in our set, for example:

- *Dickeya solani*, which will be a target species for diagnostics;
- *Pectobacterium atrosepticum*;
- *Pectobacterium carotovorum*.

It is necessary to point out that the *Dickeya* genus contains other species, which are also potato pathogens. However, the species have very little differences in their high-conservative genes and symptomatology, therefore, their distinction is not necessary.

Step 2. Selecting target nucleotide sequences.

Potential target nucleotide sequence can be found using the NCBI Nucleotide database (<https://www.ncbi.nlm.nih.gov/nuccore>), the NCBI GenBank genome annotation, or the NCBI BLAST service (<https://blast.ncbi.nlm.nih.gov>). To develop specific primer set for detection of *Dickeya* pathogens, we choose a strategy with primer design toward high-conservative housekeeping gene of histidine metabolism, called the *hisC* gene.

Two alternative types of the FASTA-header formatting are compatible with further protocol:

- `>[Taxon ID] |[Target gene ID] |[Strain ID]`
- `>[Taxon ID] |[Strain ID]`

For example, the user can insert instead of [Taxon ID] any text that will be associate with target or closely related species. Thus, for the type strain of *D. solani* is IPO2222^T the corresponding FASTA-header are the following:

- >DSO|hisC|IPO2222
- >DSO|IPO2222

The user is free to use any type of the abovementioned FASTA-header, but for high-performing screening of the target sequence might be useful a long version with gene ID. It allows to differentiate various plots and compare the data more effective. Here, we select the well-established 3-letter species abbreviation and 4-letter gene abbreviation, but any other IDs or custom text and digits can be used as well as NCBI taxon identifier or GenBank accession number. It is necessary to point out that the ID should not to contain spaces, because of the spaces cause the header gap during the alignment by Clustal Omega service.

We recommend to extract target sequence from the type strain with complete genome status. For example, the type strain of *D. solani* is IPO2222^T (GenBank accession no.: NZ_CP015137). In **Example 1.1** you can see required formatting of the FASTA-header in text (*.txt) file with nucleotide sequence of interest.

```
>DSO|hisC|IPO2222
ATGAGTATTGATAATCTGGCTCGCGACAACGTACGCCGCTGACCCCGTATCAGTCCGCTCGCCGCTGGGCGGCAACGGTGACGCTCTGGCTGAACGCCAACGAG
TACCCGAGGCGCCGAGTACCAACTGACCTTGCAGACGCTGAACCGTTACCCGGAATGTCAGCCCGTGCAGGTGATTGAACGCTATGCCGCTATGCCGGGTTG
AAACCGGAGCAGGTGCTGGTGAGCCGTGGCGCAGACGAAGGCATTGAACTGCTGATTTCGCGCGTTCTGCGAACCAGGGTAAAGACGCCATTCTGTCTGTCCACCG
ACCTACGGCATGTACGCCGTAGCGCTGAAACGTTTCGGCGTAGAGCGCCGCTAGTGCCTCGACCGCGACTGGCAACTGGATCTGGCCGCGATCGAAGGCGCG
CTGGACAACGTTAAAGTGATTATGTTTGTTCGCCCAACAACCCGACCGGAAACCTGATCAACCCAGAGGATCTGCGCCGCTCTGCTGGAGTTGGCGCGCGGTCGG
GCCATCGTCGCTATCGACGAAGCCTATATTGAGTTTTGCCCGCAGGCGACACGGTCGGCTGGCTGGCGGAATACCCGCATCAGGTCGTATTGCGCACGCTGTCC
AAAGCCTTCGCGCTGGCCGGCCTGCGTTGTGGTTTCACCCCTGCCAGCGCAGGGGTGATTCAACTGTGCTCAAAGTGATCGCCCCCTACCCACTGGCGTTGCCG
GTAGCGGATATCGCCGCGCAGGCGCTGAGCGCTGAGGGTCTGGCGGACATGCGCCGCAATGTGGACGAGGTGCGGGAAACCCCGCTGGCTGCGTGAGTCGCTG
AAAGCGTTGTCGAGCATAGAAACGGTATACGCCAGCGAAAGCAATTATCTACTGGTGCGCTTCGCCGATTTCGCCAACCGTCTTCAAAACCTTGTGGGATCAGGGC
ATTATTTTACGCGATCAGAACAAACAGCCGGGGCTGGCGGGTTGTCTGCGCATTACCATCGGCAACCGCTACGAGTGCGAGCGCGTCATCAGCGCATTGCAGGCG
TTATCCGGCAAAACGGCATAA
```

Example 1.1. Required formatting of the target sequence with modified FASTA-header.

Step 3. Exporting orthologous sequences.

Using the NCBI BLAST service and the target nucleotide sequence as a reference, the user should to export some orthologous sequences from closely related bacterial species. In present example, the *hisC* sequences were extracted from three bacterial species, each of them was represented by five different strains, including type strain. All sequences with modified FASTA-header should be merged in one text (*.txt) file as shown below in **Example 1.2**.

```

>DSO|hisC|IPO2222
ATGAGTATTGATAAATCTGGCTCGCGACAACGTACGCCGCTGACCCCGTATCAGTCCGCTCGCCGCTGGGCGGCAACGGTGACGCTCTGGCTGAACGCCAACGAG
TACCCGCGAGGCGCGCAGTACCAACTGACCTTGACAGCGCTGAACCGTTACCCGGAATGTACAGCCCGTGCAGGTGATTGAACGCTATGCCGCCATGCCGGGTTG
AAACCGGAGCAGGTGCTGGTGAAGCGTGGCGCAGACGAAGGCATTGAACTGCTGATTGCGCGCTTCTGCGAACCAGGGTAAAGACGCCATTCTGTTCTGTCCACCG
ACCTACGGCATGTACGCCGTGAGCGCTGAAACGTTTCGGCGTAGAGCGCCGCTAGTGCCCTCGACCGCCGACTGGCAACTGGATCTGGCCGCGATCGAAGGCGCG
CTGGACAACGTAAAGTGATTATGTTTGTTCGCCCAACAACCCGACCGGAAACCTGATCAACCCAGAGGATCTGCGCGCTCTGCTGGAGTTGGCGCGCGGTGCG
GCCATCGTCGCTATCGACGAAGCCTATATTGAGTTTTCGCCGACGCGACACGGTTCGGCTGGCTGGCGGAATACCCGCATCAGGTGCTATTGCGCACGCTGTCC
AAAGCCTTCGCGCTGCGCGGCTGCGTTGTGGTTTACCCTTGCCAGCGCAGGGGTGATTCAACTGCTGCTCAAAGTGATCGCCCCCTACCCACTGGCGTTGCCG
GTAGCGGATATCGCCGCGCAGGCGCTGAGCGCTGAGGGTCTGGCGGACATGCGCCGCAATGTGGACGAGGTGCGGGAAAACCGCCGCTGGCTGCGTGAGTCGCTG
AAAGCGTTGTGCGAGCATAGAACCGGTATACGCCAGCGAAAGCAATTATCTACTGGTGCCTTCGCCGATTTCGCCAACCGTCTTCAAAACCTTGTGGGATCAGGGC
ATTATTTTACGCGATCAGAACAACAGCCGGGGCTGGCGGGTTGTCTGCGCATTACCATCGGCAACCGCTACGAGTGCGAGCGCGTCATCAGCGCATTCAGGGC
TTATCCGGCAAAACGGCATAA

>DSO|hisC|F012
ATGAGTATTGATAAATCTGGCTCGCGACAACGTACGCCGCTGACCCCGTATCAGTCCGCTCGCCGCTGGGCGGCAACGGTGACGCTCTGGCTGAACGCCAACGAG
TACCCGCGAGGCGCGCAGTACCAACTGACCTTGACAGCGCTGAACCGTTACCCGGAATGTACAGCCCGTGCAGGTGATTGAACGCTATGCCGCCATGCCGGGTTG
AAACCGGAGCAGGTGCTGGTGAAGCGTGGCGCAGACGAAGGCATTGAACTGCTGATTGCGCGCTTCTGCGAACCAGGGTAAAGACGCCATTCTGTTCTGTCCACCG
ACCTACGGCATGTACGCCGTGAGCGCTGAAACGTTTCGGCGTAGAGCGCCGCTAGTGCCCTCGACCGCCGACTGGCAACTGGATCTGGCCGCGATCGAAGGCGCG
CTGGACAACGTAAAGTGATTATGTTTGTTCGCCCAACAACCCGACCGGAAACCTGATCAACCCAGAGGATCTGCGCGCTCTGCTGGAGTTGGCGCGCGGTGCG
GCCATCGTCGCTATCGACGAAGCCTATATTGAGTTTTCGCCGACGCGACACGGTTCGGCTGGCTGGCGGAATACCCGCATCAGGTGCTATTGCGCACGCTGTCC
AAAGCCTTCGCGCTGCGCGGCTGCGTTGTGGTTTACCCTTGCCAGCGCAGGGGTGATTCAACTGCTGCTCAAAGTGATCGCCCCCTACCCACTGGCGTTGCCG
GTAGCGGATATCGCCGCGCAGGCGCTGAGCGCTGAGGGTCTGGCGGACATGCGCCGCAATGTGGACGAGGTGCGGGAAAACCGCCGCTGGCTGCGTGAGTCGCTG
AAAGCGTTGTGCGAGCATAGAACCGGTATACGCCAGCGAAAGCAATTATCTACTGGTGCCTTCGCCGATTTCGCCAACCGTCTTCAAAACCTTGTGGGATCAGGGC
ATTATTTTACGCGATCAGAACAACAGCCGGGGCTGGCGGGTTGTCTGCGCATTACCATCGGCAACCGCTACGAGTGCGAGCGCGTCATCAGCGCATTCAGGGC
TTATCCGGCAAAACGGCATAA

...

>PCA|hisC|YCD57
ATGAGCAGCATTGAAGAACTGGCACGTGCCAACGTCCGTGCGTTGACCCCTACCAATCTGCCCGTCGTCTGGGCGGTAACGGCGATGTCTGGTTGAATGCGAAT
GAATATCCTGAAGCACCAGAGTTTCAGTTAACGCTGCAACCGCTGAACCGCTACCCGGAGTGTCAGCCGGTACCGGTGATCAACCGTTACGCTGAATATGCTGGC
GTAACGCCGAGCAGGTTCCTGGTCAGCCGCGGTGCCGACGAAGGTATTGAGCTGCTGATTGCGCGCTTCTGCGAACCAGGGGAAAGATGCCATCCTGTTCTGTCCG
CCAACCTACGGCATGTACGCCGTGAGCGCTGAGACGTTTGGCGTCGAACGTCGCACCGTAGCCAGTAAATCAGACTGGCAGTTGGATCTCGATGCGATTGAAGCA
CAACTGGACGGCACCAGAGTATTTACGTATGCAGCCAAATAACCCAACCGGCAACCTGATTGCACGGGAGGATTTACGCCAGTTGCTCACAAATGGCACAGGGC
AAAGCGTTGGTCGTTATCGATGAGGCTTACATCGAGTTCTGTCCGAGGCGTCCACCGCAGTCTGGCTAAGCGAGTTTCTCATTTGGTGATTCTGCGTACCTTG
TCCAAAGCCTTTTCACTCGCGGGGCTGCGCTGTGGCTTACGCTGGCGAACCTGAAGTGATTACAGTTCTGCTGAAGGTGATTGCCCTTATCCGCTGTCTACT
CCGGTGGCAGATATCGCCGCGCAGGCGCTGAGTCACGAAGGATCGCCAAATGAAGGCGAATGTTGAGGAAATTACCTCCGCTCGCCGCTGGCTGAGCGATGCC
TTAAAGACATTCCCTGCGTTGAAGAAATCTTCCCGAGCGAGGCAACTATTGCTGGTGCCTTTACCGCTTCCCTTCTGTATTTAAACGCTGTGGGATCAA
GGCATTATTCTGCGTGACCAAAATAAGCAACCGAGCCTCGCAGGCTGTTTGGCGATTACTATCGGCAACCGCTACGAATGTGAACGCGTTGTTGCTGCATTACAA
TCATTGCGGGGCATCAACGCTTAA

```

Example 1.2. Single FASTA-file of orthologous sequences with modified headers.

Step 4. Making multiple sequence alignment.

To perform multiple alignment of the orthologs, we recommend to use the Clustal Omega service (<https://www.ebi.ac.uk/Tools/msa/clustalo>) with parameters that shown in **Figure 1.1**. It is necessary to select an "Input" value of the "Order" menu to save the species order in single FASTA-file. This point is important because the MorphoCatcher algorithm is based to analyze the strains of different species as separate groups and find their conservative species-specific mutations. Then the user should to save an alignment text (*.txt) file as shown in **Example 1.3**.

EMBL-EBI

Services

Research

Training

Industry

About us

Q

EMBL-EBI

Hinxton

Clustal Omega

Input form

Web services

Help & Documentation

Bioinformatics Tools FAQ

Feedback

Share

STEP 1 - Enter your input sequences

Enter or paste a set of

DNA

sequences in any supported format:

>DSO|hisC|IPO2222

ATGAGTATTGATAATCTGGCTCGCGACAACTACGCCGCTGACCCCGTATCAGTCCGCTCGCCGCTGGGCGGCAACGGTGACGTCTGGCTGAACGCCA

ACGAGTACCCGCGAGGCGCCGAGTACCACTGACCTTGACAGCGCTGAACCGTTACCCGGAATGTCAGCCCGTGCAGGTGATTGAACGCTATGCCGCCTA

TGCCGGGTTGAAACCGGAGCAGGTGCTGGTGAGCCGTGGCGCAGACGAAGGCATTGAAGTCTGATTTCGCGCGTTCTGCGAACCAGGTAAGACGCCAT

TCTGTTCTGTCCACCGACCTACGGCATGTACGCCGTGAGCGCTGAAACGTTGCGCGTAGAGCGCCGCGTAGTGCCCTCGACCGCCGACTGGCAACTGGA

TCTGGCCGCGATCGAAGGCGCGCTGGACAACGTTAAAGTGATTATGTTTGTTCGCCCAACAACCCGACCGGAAACCTGATCAACCCAGAGGATCTGCGC

CGTCTGCTGGAGTTGGCGCGCGGTGCGGCCATCGTCGCTATCGACGAAGCCTATATTGAGTTTGGCCGCGGCGACACCGGTGCGCTGGCTGGCGGAA

TACCCGCATCAGGTGCTATTGCGCACGCTGTCCAAAGCCTTGCAGCTGGCCGCGCTGCGTTGTGGTTTACCCTTGCCAGCGCAGGGGTGATTCAACTGC

Or, upload a file:

Обзор...

Файл не выбран.

[Use a example sequence](#) | [Clear sequence](#) | [See more example inputs](#)

STEP 2 - Set your parameters

OUTPUT FORMAT

ClustalW with character counts

DEALIGN INPUT SEQUENCES	MBED-LIKE CLUSTERING GUIDE-TREE	MBED-LIKE CLUSTERING ITERATION	NUMBER of COMBINED ITERATIONS
no	yes	yes	default(0)
MAX GUIDE TREE ITERATIONS	MAX HMM ITERATIONS	ORDER	
default	default	input	

STEP 3 - Submit your job

☐ Be notified by email *(Tick this box if you want to be notified by email when the results are available)*

Submit

Figure 1.1. Screenshot of the Clustal Omega service with required options for alignment. The "Input" option of "Order" menu should be selected.

CLUSTAL O(1.2.4) multiple sequence alignment

```

DSO|hisC|IPO2222      ---ATGAGTATTGATAAATCTGGCTCGCGACAACGTACGCCGCCTGACCCCGTATCAGTCC      57
DSO|hisC|F012         ---ATGAGTATTGATAAATCTGGCTCGCGACAACGTACGCCGCCTGACCCCGTATCAGTCC      57
DSO|hisC|MK10         ---ATGAGTATTGATAAATCTGGCTCGCGACAACGTACGCCGCCTGACCCCGTATCAGTCC      57
DSO|hisC|MK16         ---ATGAGTATTGATAAATCTGGCTCGCGACAACGTACGCCGCCTGACCCCGTATCAGTCC      57
DSO|hisC|RNS08.23.3.1.A ---ATGAGTATTGATAAATCTGGCTCGCGACAACGTACGCCGCCTGACCCCGTATCAGTCC      57
PAT|hisC|NCPB549      ATGAGCAGCATTGAAGAAGTGGCAGCTGCCAACGTCCGTGCATTGACACCTTATCAGTCT      60
PAT|hisC|JG10-08      ATGAGCAGCATTGAAGAAGTGGCAGCTGCCAACGTCCGTGCATTGACACCTTATCAGTCT      60
PAT|hisC|NCPB3404     ATGAGCAGCATTGAAGAAGTGGCAGCTGCCAACGTCCGTGCATTGACACCTTATCAGTCT      60
PAT|hisC|SCRI1043     ATGAGCAGCATTGAAGAAGTAGCACGCGCCAACGTCCGTGCATTGACACCGTATCAGTCT      60
PAT|hisC|21A          ATGAGCAGCATTGAAGAAGTGGCAGCTGCCAACGTCCGTGCATTGACACCTTATCAGTCT      60
PCA|hisC|NCPB312      ATGAGCAGCATTGAAGAAGTGGCAGCTGCCAACGTCCGTGCATTGACACCGTATCAGTCT      60
PCA|hisC|BCT2         ATGAGCAGTATTGAAGAAGTGGCAGCGCCAATGTCCGTGCGCTGACGCGGTACCAATCT      60
PCA|hisC|BCT5         ATGAGCAGCATTGAAGAAGTGGCAGCGCCAATGTCCGTGCGCTGACGCGGTACCAATCT      60
PCA|hisC|WPP14        ATGAGCAGCATTGAAGAAGTGGCAGCGCCAATGTCCGTGCGCTGACGCGGTACCAATCT      60
PCA|hisC|YCD57        ATGAGCAGCATTGAAGAAGTGGCAGCTGCCAACGTCCGTGCGTTGACCCCTACCAATCT      60
                        *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *
...

DSO|hisC|IPO2222      GATCAGAAACAAACAGCCGGGGTGGCGGGTTGTCTGCGCATTACCATCGGCAACCGCTAC      1017
DSO|hisC|F012         GATCAGAAACAAACAGCCGGGGTGGCGGGTTGTCTGCGCATTACCATCGGCAACCGCTAC      1017
DSO|hisC|MK10         GATCAGAAACAAACAGCCGGGGTGGCGGGTTGTCTGCGCATTACCATCGGCAACCGCTAC      1017
DSO|hisC|MK16         GATCAGAAACAAACAGCCGGGGTGGCGGGTTGTCTGCGCATTACCATCGGCAACCGCTAC      1017
DSO|hisC|RNS08.23.3.1.A GATCAGAAACAAACAGCCGGGGTGGCGGGTTGTCTGCGCATTACCATCGGCAACCGCTAC      1017
PAT|hisC|NCPB549      GACCAAAATAAGCAACCCAGCCTCGCAGGCTGTTTGCGCATTACTATCGGCAACCGCTAC      1020
PAT|hisC|JG10-08      GACCAAAATAAGCAACCCAGCCTCGCAGGCTGTTTGCGCATTACTATCGGCAACCGCTAC      1020
PAT|hisC|NCPB3404     GACCAAAATAAGCAACCCAGCCTCGCAGGCTGTTTGCGCATTACTATCGGCAACCGCTAC      1020
PAT|hisC|SCRI1043     GACCAAAATAAGCAACCCAGCCTCGCAGGCTGTTTGCGCATTACTATCGGCAACCGCTAC      1020
PAT|hisC|21A          GACCAAAATAAGCAACCCAGCCTCGCAGGCTGTTTGCGCATTACTATCGGCAACCGCTAC      1020
PCA|hisC|NCPB312      GACCAAAATAAGCAACCCAGCCTCGCAGGCTGTTTGCGCATTACTATCGGCAACCGCTAC      1020
PCA|hisC|BCT2         GACCAAAATAAGCAACCCAGCCTCGCAGGCTGTTTGCGCATTACTATCGGCAACCGCTAC      1020
PCA|hisC|BCT5         GACCAAAATAAGCAACCCAGCCTCGCAGGCTGTTTGCGCATTACTATCGGCAACCGCTAC      1020
PCA|hisC|WPP14        GACCAAAATAAGCAACCCAGCCTCGCAGGCTGTTTGCGCATTACTATCGGCAACCGCTAC      1019
PCA|hisC|YCD57        GACCAAAATAAGCAACCCAGCCTCGCAGGCTGTTTGCGCATTACTATCGGCAACCGCTAC      1020
                        *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *

DSO|hisC|IPO2222      GAGTGCGAGCGCGTCATCAGCGCATTGCAGGCGTTATCCGGCAAAACGGCATAA      1071
DSO|hisC|F012         GAGTGCGAGCGCGTCATCAGCGCATTGCAGGCGTTATCCGGCAAAACGGCATAA      1071
DSO|hisC|MK10         GAGTGCGAGCGCGTCATCAGCGCATTGCAGGCGTTATCCGGCAAAACGGCATAA      1071
DSO|hisC|MK16         GAGTGCGAGCGCGTCATCAGCGCATTGCAGGCGTTATCCGGCAAAACGGCATAA      1071
DSO|hisC|RNS08.23.3.1.A GAGTGCGAGCGCGTCATCAGCGCATTGCAGGCGTTATCCGGCAAAACGGCATAA      1071
PAT|hisC|NCPB549      GAATGTGAACGCGTTGTTGCTGCATTACAATCATTTGCCGGGCATCAACGCTTAA      1074
PAT|hisC|JG10-08      GAATGTGAACGCGTTGTTGCTGCATTACAATCATTTGCCGGGCATCAACGCTTAA      1074
PAT|hisC|NCPB3404     GAATGTGAACGCGTTGTTGCTGCATTACAATCATTTGCCGGGCATCAACGCTTAA      1074
PAT|hisC|SCRI1043     GAATGTGAACGCGTTGTTGCTGCATTACAATCATTTGCCGGGCATCAACGCTTAA      1074
PAT|hisC|21A          GAATGTGAACGCGTTGTTGCTGCATTACAATCATTTGCCGGGCATCAACGCTTAA      1074
PCA|hisC|NCPB312      GAATGTGAACGCGTTGTTGCTGCATTACAATCATTTGCCGGGCATCAACGCTTAA      1074
PCA|hisC|BCT2         GAATGTGAACGCGTTGTTGCTGCATTACAATCATTTGCCGGGCATCAACGCTTAA      1074
PCA|hisC|BCT5         GAATGTGAACGCGTTGTTGCTGCATTACAATCATTTGCCGGGCATCAACGCTTAA      1074
PCA|hisC|WPP14        GAATGTGAACGCGTTGTTGCTGCATTACAATCATTTGCCGGGCATCAACGCTTAA      1073
PCA|hisC|YCD57        GAATGTGAACGCGTTGTTGCTGCATTACAATCATTTGCCGGGCATCAACGCTTAA      1074
                        *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *  *

```

Example 1.3. Multiple sequence alignment of the *hisC* orthologous sequences. Conservative nucleotide positions are highlighted by "*" symbols.

Step 5. Masking species-specific mutations.

Here, we suggest to use the MorphoCatcher service (<http://morphocatcher.ru>) for processing of the MSA. The MorphoCatcher allows the user to upload the alignment and select the output file format, which depends on computer operating system (e.g. Windows or UNIX). The default parameter set is "Windows" output file format.

After alignment uploading the user should to click on the "Catch" button, which starts processing of data to detect species-specific mutations (e.g. single nucleotide polymorphisms, insertions, and deletions). Then, the new tab with results will shows a graphic representation of the mutation density using a sliding-window function (**Figure 1.2**).




Figure 1.2. Screenshot of the MorphoCatcher service with option of the output format.

RESULTS

ID18022019224247




Figure 1.3. Screenshot of tab with results after processing of the multiple sequence alignment.

Among output files (**Figure 1.3**) the user can find the following files with useful data:

- Processed MSA with masked species-specific mutations;
- Sliding-window plot of gene polymorphisms for each species (**Figure 1.4**);
- Ready-to-use files for specific primer design using the PrimerExplorer service.

All these files can be download as ZIP-archive with corresponding identification number of the project. Graphic files have a high resolution for publications and available with different extensions (*.eps, *.pdf, *.png, *.svg, and *.tiff). Some large files (*.tiff, *.eps, and *.svg) are included only in ZIP-archive. Other graphic (*.pdf and *.png) and text (*.txt) files are available for download and preview in new tab of your browser.




Figure 1.4. Sliding-window plot of the *hisC* gene polymorphism in closely related species. Peaks correspond to the high density of species-specific mutations.

The sequence regions with maximal average mutation index correspond to the high density of species-specific mutations. To select most suitable target region, we recommend the average mutation index values from 2 and higher. If the target selected will be compatible with the design algorithm of PrimerExplorer, the most of species-specific mutations can be included at the primer ends for better specificity.

The user can also inspect positions of the species-specific mutations manually by using processed alignment text file from the output. For this procedure the processed alignment file contains number of species-specific mutations for each 20 nucleotides (**Example 1.4**).

UPLOAD FILE: DSO.hisC.IPO2222.txt

```

1  ATGAGTATTG ATAATCTGGC TCGCGACAACT GTACGCCGCG TGACCCCGTA TCAGTCGCGT CGCCGCGCTGG GCGGCAACGG  80
   *-----*
81  TGACGTCTGG CTGAACGCCA ACGAGTACCC GCAGCGCGCG CAGTACCAAC TGACCTTGCA GACGCTGAAC CGTTACCCGG  160
   *-----*
161 AATGTCAGCC CGTGCAGGTG ATTGAACGCT ATGCCGCCTA TGCCGGGTTG AACCGGAGC AGGTGCTGCT GAGCCGTGGC  240
   *-----*
241 GCAGACGAAG GCATTGAAC TCTGATTGCG GCGTTCTGCG AACCGGGTAA AGACGCCATT CTGTTCTGTC CACCGACCTA  320
   *-----*
321 CGGCATGTAC GCCGTCAGCG CTGAAACGTT CGGCGTAGAG CGCCGCGTAG TGCCCTCGAC CGCCGACTGG CAACTGGATC  400
   *-----*
401 TGGCCGCGAT CGAAGGCGCG CTGGACAACG TTAAAGTGAT TTATGTTTGT TCGCCCAACA ACCCGACCGG AAACCTGATC  480
   *-----*
481 AACCCAGAGG ATCTGCGCGC TCTGCTGGAG TTGGCGCGCG GTCGGGCGAT CGTCGCTATC GACGAAGCCT ATATTGAGTT  560
   *-----*
561 TTGCCCGCAG GCGACCAACG TCGGCTGGCT GCGCGAATAC CCGCATCAGG TCGTATTGCG CACGCTGTCC AAAGCCTTCG  640
   *-----*
641 CGCTGGCCGG CCTGCGTTGT GGTTCACCC TTGCCAGCGC AGGGGTGATT CAATGCTGCT TCAAAGTGAT CGCCCCCTAC  720
   *-----*
721 CCACTGGCGT TGCCGGTAGC GGATATCGCC GCGCAGGCGC TGAGCGCTGA GGGTCTGGCG GACATGCGCC GCAATGTGGA  800
   *-----*
801 CGAGGTGCGG GAAACCGGCC GCTGGCTGCG TGAGTCGCTG AAAGCGTTGT CGAGCATAGA AACGGTATAC GCCAGCGAAA  880
   *-----*

```

Number of Primer Candidates: F1=194, F2=133, F3=245, B1=164, B2=140, B3=260, FIP=228, BIP=204
812 Primer set(s) were generated.

Set Mutation
Mut/Cons
Clear

Fixed Primer
F3
F2
F1
B1
B2
B3
Clear

Save Target

Design Option
☐ Default
☐ Common
☒ Specific

1. Select Range
☒ Ignore range
☐ Within F2-B2
☐ Between F1c-B1c Targeting Range

2. Generate
Generate 812 sets were generated.

3. Display
Display Page 1 Displayed. Sorting Rule Mutation(s)[descending order]

If you can move to "Basic Designing", please click below.

Basic Designing

Parameter Condition Normal Save Parameter Reset Parameter

Figure 1.5. Interface of the PrimerExplorer service after uploading target nucleotide sequence. Any specific mutation that can be covered by primer ends is highlighted by "-" symbol.

To design primer sets with species-specific mutations at 5'- or 3'-ends, the user should to configure these important primer parameters of the PrimerExplorer service:

- To select "Specific" mode of the "Primer option" menu;
- To click on the "Detail settings" button to turn on an expert mode;
- To select a required GC-content using the "Parameter condition" menu;
- To click on the "Generate" button;
- To set "Descending order" of primer sets in the "Sorting rule" menu;
- To click on the "Display" button.

After that the user will be redirected to the page with various combinations of primer annealing sites, which can cover various species-specific mutations (**Figure 1.6**). Some of the mutations are localized at primer ends (**Figure 1.7**) that able to increase specificity of the amplification.




Figure 1.6. The PrimerExplorer output page with full view of designed primer sets. Detected species-specific mutations for strains of *D. solani* are highlighted by red color. The same mutations can be covered by oligos from different primer sets.




Figure 1.7. Detailed view of the primer set sequences with different thermodynamic parameters. Species-specific mutations that covered by primer ends are highlighted by red color.

Case study 2: the viral targets and whole-genome multiple alignment

Step 1. Defining strains from different serotypes.

Well-known viral potato pathogen is potato virus Y (PVY) that causes mosaic and leaf drop. Genome of PVY is the single-strand (+) RNA molecule and the length is approximately 9.7 kb. Most potato infections are caused by two main strains that have similar genome organization:

- PVY^O — the ordinary strain;
- PVY^N — the necrotic strain.

There are also various new recombinants derived from the abovementioned strains [2]. All these strains and recombinants usually belong to two serotypes, called O- and N-serotype. Due to different symptoms caused by the PVY serotypes, it is necessary to distinguish them for the potato seed certification purposes. In this case study, we attempt to use whole genome sequences for designing LAMP primer sets for the detection of different PVY serotypes. The genome size is applicable for fast multiple alignment of different PVY strains. Our main goal was to decide which genome region can be used for the developing of serotype-specific assays.

Step 2. Selecting target genome sequences.

Here, we attempt to distinguish the PVY serotypes using the MSA of original PVY strains and their recombinants and try to apply the strategy of serotype-specific polymorphisms search.

Step 3. Exporting genome sequences.

To obtain whole-genome sequences of PVY strains, we collect data from NCBI Nucleotide database (<https://www.ncbi.nlm.nih.gov/nuccore>) and scientific literature [3] by entering corresponding accession numbers in the service search bar. The PVY genomes with corresponding accession numbers can be easily accessed at the link:

<https://www.ncbi.nlm.nih.gov/nuccore/AJ585195,AJ585196,AJ890349,DQ157178,EF558545,HQ912864,HQ912865,HQ912888,U09509,X12456,AJ585197,AJ585342,AJ889866,AY166866,AY884982,AY884983,FJ666337,HQ912869,M95491,X97895>

The user can to download all these genome sequences as single FASTA-file using corresponding menu of NCBI Nucleotide web page (see **Figure 2.1**) and after that to change the *.fasta extension to *.txt by renaming of the file for comfort data processing. For correct protocol execution the user should to modify the FASTA-header for each genome of interest. The user is free to include any additional service information to the FASTA-header. We use the accession number, because of the number is already part of the standard FASTA-header of each genome at the NCBI Nucleotide database. Thus, we have two main groups that correspond to O- and N-serotype, therefore, the FASTA-headers for each PVY genome can be modified using the following mask:

- >PVYO| [Accession]
- >PVYN| [Accession]

The final view of PVY genome sequences is shown in **Example 2.1**.

Items: 20

- ☐ [Potato virus Y mRNA for polyprotein \(pol gene\), isolate SASA-110](#)
1. 9,651 bp linear RNA
Accession: AJ585195.1 GI: 53913350
[Protein](#) [Taxonomy](#)
[GenBank](#) [FASTA](#) [Graphics](#)
- ☐ [Potato virus Y mRNA for polyprotein \(pol gene\), isolate SCRI-O](#)
2. 9,698 bp linear RNA
Accession: AJ585196.1 GI: 53913352
[Protein](#) [Taxonomy](#)
[GenBank](#) [FASTA](#) [Graphics](#)
- ☐ [Potato virus Y strain O gene for polyprotein, genomic RNA, isolate LW](#)
3. 9,699 bp linear RNA
Accession: AJ890349.1 GI: 90968473
[Protein](#) [PubMed](#) [Taxonomy](#)
[GenBank](#) [FASTA](#) [Graphics](#)

☒ Complete Record
☐ Coding Sequences
☐ Gene Features

Choose Destination
☒ File ☐ Clipboard
☐ Collections

Download 20 items.

Format
FASTA ▾

Sort by
Accession ▾

Show GI ☐

Create File

Recent activity [Turn Off](#) [Clear](#)

Figure 2.1. The NCBI Nucleotide web page with list of well-known PVY genomes. To download all sequences, the user should click the "Send to:" menu, then to select "Complete record" and "File" destination, and choose the "FASTA" format with sorting by accession number. To save a text file, click on "Create file" button to save a text file.

```
>PVYO|AJ585195
AAACCAACGCAAAAACACTTATAAACGCTTATCTCACTCAAGCAACTTACTAAGTTTCAGTTTAAATCA
TTTCCTTGCAATCTCTTAAACAATATTGGAACACTATTTCAACTCAACAAGCAATTTTCATCACTTCCAAC
CAATTTTAGATCCTCAATGGCAACTTACATGTCAACAATCTGTTTCGGTTCGTTGAATGCAAGCTACCA
TACTCACCCGCTGTTGCGGGCATATTGTGAAGGAACGAGAAGTGCTGGCTTCCGTTGATCCTTTTCGCAG
ATCTGGAACACAACTTAGTGACGATTGCTCAAGCAAAAATATGCTACTGTTGCTGTGCTCAAGAACGG
TACTCTTACGTACCGATACAAGAATGATGCCAGATAACGCGCATCCAGAAGAACTGGAAAGGAAGGAT
AGGGAAGAATATCACTTCCAGATGGCAGCTCCTAGTATTGTGTCAAAAATTACTATAGCTGGTGGAGATC
CTCCATCAAAGTCAGAGTCACAAGCACAAGAGGTGTCATTATACAACTCCAAGGGTGCGTAAAGTCAA
GACACGCCCCATAACAAAGTTGACAGAAGGCCAGATGAATCATCTCATTAAGCAGGTGAAGCAGATTATG
TCGGAGAAGAGAGGGTCCGTCCACTTAATTAGTAAGAAGACCCTCATGTTCAATATAAGGAGATACTTG
...
TGAGTACCCGTTGAAACCAATTGTTGAGAATGCAAAACCAACCCCTTAGGCAAAATCATGGCACATTTCTCA
GATGTTGAGAGCGTATATAGAAATGCGCAACAAAAGGAACCATATATGCCACGATATGGTTTAATTC
GAAATCTGCGGGATGTGGGTTTAGCGCGCTATGCCTTTGACTTCTACGAGGTACATCACGAACACCACT
GAGGGCTAGGGAAGCGCACATTCAAATGAAGCGCGCAGCATTTGAAATCAGCCCAACCTCGACTTTTCGGG
TTGGACGGTGGCATCAGTACACAAGAGGAGAACACAGAGAGGCACACCACCGAGGATGTCCTCCAAAGTA
TGCATACTCTACTTGGAGTCAAGAACATGTGATGTAGTGTCTCTCCGACGATATATAAGTATTACATA
TGCAGTAAGTATTTTGCTTTTCTGTACTACTTTTATCATAATTAATAATCAGTTTGAATATTACTAAT
AGATAGAGGTGGCAGGGTGATTTCGTGATTGTGGTGACTCTATCTGTTAATTTTCGATTATTAAGTCTTA
GATAAAGTGCCGGGTTGTCGTTGTTGGATGATTATCGATTAGGTGATGTTGCGATTCTGTCGTAGC
AGTGACTATGTCTGAATCTATCTGCTTGGGTGGTGTGATTTCGTCATGACAGTGACTG
```

Example 2.1. The single FASTA-file of PVY genome sequences with modified headers.

Step 4. Making multiple sequence alignment.

Due to the PVY genome size, it is possible to obtain the MSA of PVY genomes by using the Clustal Omega service (<https://www.ebi.ac.uk/Tools/msa/clustalo>) with the same parameters that shown in **Figure 1.1**.

CLUSTAL O (1.2.4) multiple sequence alignment

```

PVYO|AJ585195 -----AAACCAACGCAAAAAACACTTATAAACGCTTA 31
PVYO|AJ585196 -AATTAACAACACTCAATACAACATAAGAAAAACAACGCAAAAAACACTCATAAACGCTTA 59
PVYO|AJ890349 AAATTAACAACACTCAATACAACATAAGAAAAACAACGCAAAAAACACTCATAAACGCT-- 58
PVYO|DQ157178 -----CAACGCAAAAAACACTCATAAACGCTTT 27
PVYO|EF558545 -AATTAACAACACTCAATACAACATAAGAAAAACAACGCAAAAAACACTCATAAACGCTCA 59
PVYO|HQ912864 -AATTAACAACACTCAATACAACATAAGAAAAACAACGCAAAAAACACTCATAAACGCTTA 59
PVYO|HQ912865 -AATTAACAACACTCAATACAACATAAGAAAAACAACGCAAAAAACACTCATAAACGCTTA 59
PVYO|HQ912888 -----TCAACGCAAAAAACACTCATAAAAGCTTA 28
PVYO|U09509 -AATTAACAACACTCAATACAACATAAGAAAAACAGCGCAAAAAACACTCATAAACGCTTA 59
PVYO|X12456 -AATTAACAACACTCAATACAACATAAGAAAAACAACGCAAAAAACACTCATAAACGCTCA 59
PVYN|AJ585197 -----AAATCAACGCAAAAAACACTCACAAAAGCTTT 31
PVYN|AJ585342 -----GCAAAAAACACTCACAAAAGCTTT 23
PVYN|AJ889866 AAATTAACAACACTCAATACAACATAAGAAAAACAACGCAAAAAACACTCATAAACGCTTA 60
PVYN|AY166866 -AATTAACAACACTCAATACAACATAAGAAAAACAACGCAAAAAACACTCACAAAAGCTTT 59
PVYN|AY884982 -AATTAACAACACTCAATACAACATAAGAAAAACAACGCAAAAAACACTCACAAAAGCTTT 59
PVYN|AY884983 -AATTAACAACACTCAATACAACATAAGAAAAACAACGCAAAAAACACTCACAAAAGCTTT 59
PVYN|FJ666337 ----- 0
PVYN|HQ912869 -----TCAACGCAAAAAACACTCATAAAAGCTTT 28
PVYN|M95491 AAATTAACAACACTCAATACAACATAAGAAAAACAACGCAAAAAACACTCACAAAAGCTTT 60
PVYN|X97895 AAATTAACAACACTCAATACAACATAAGAAAAACAACGCAAAAAACACTCACAAAAGCTTT 60

...

PVYO|AJ585195 ATCGATTAGGTGATGTTGCGATTCTGTCTAGCAGTGACTATGTCTGAATCTATCTGCTT 9617
PVYO|AJ585196 ATCGATTAGGTGATGTTGCGATTCTGTCTAGCAGTGACTATGTCTGGATCTATCTACTT 9645
PVYO|AJ890349 ATCGATTAGGTGATGTTGCGATTCTGTCTAGCAGTGACTATGTCTGGATCTATCTACTT 9646
PVYO|DQ157178 ATCGATTAGGTGGTGTGCGATTCTGTCTAGCAGTGACTATGTCTGGATCCATCTGCTT 9617
PVYO|EF558545 ATCGATTAGGTGATGTTGCGATTCTGTCTAGCAGTGACTATGTCTGGATCTATCTACTT 9645
PVYO|HQ912864 ATCGATTAGGTGATGTCGCGATTCTGTCTAGCAGTGACTATGTCTGGATCTATCTGCTT 9645
PVYO|HQ912865 ATCGATTAGGTGATGTTGCGATTCTGTCTAGCAGTGACTATGTCTGGATCTATCTGCTT 9645
PVYO|HQ912888 ATCGATTAGGTGATGTTGCGATTCTGTCTAGCAGTGACTATGTCTGGATCTATCTGCTT 9614
PVYO|U09509 ATCGATTAGGTGATGTTGCGATTCTGTCTAGCAGTGACTATGTCTGGATCTATCTGCTT 9645
PVYO|X12456 ATCGATTAGGTGATGCTGTGATTCTGTCTAGCAGTGACTATGTCTGGATTTAGTTACTT 9651
PVYN|AJ585197 ATCGATTAGGTGATGTTGCGA-TTGTCTAGCAGTGACTATGTCTGGATTTAGTTACTT 9619
PVYN|AJ585342 ATCGATTAGGTGATGTTGCGATTCTGTCTAGCAGTGACTATGTCTGGATCTATCTGCTT 9613
PVYN|AJ889866 ATCGATTAGGTGATGTTGCGATTCTGTCTAGCAGTGACTATGTCTGGATCTATCTACTT 9646
PVYN|AY166866 ATCGATTGGGTGATGTTGCGA-TTGTCTAGCAGTGACCATGTCTGGATTTAGTTACTT 9647
PVYN|AY884982 ATCGATTAGGTGATGTTGCGATTCTGTCTAGCAGTGACTATGTCTGGATCTATCTGCTT 9649
PVYN|AY884983 ATCGATTAGGTGATGTTGCGA-TTGTCTAGCAGTGACTATGTCTGGATTTAGTTACTT 9647
PVYN|FJ666337 ATCGATTAGGTGATGTTGCGATTCTGTCTAGCAGTGACTATGTCTGGATCTATCTACTT 9558
PVYN|HQ912869 ATCGATTAGGTGATGTTGCGATTCTGTCTAGCAGTGACTATGTCTGGATCTATCTGCTT 9618
PVYN|M95491 ATCGATTAGGTGATGTTGCGATTCTGTCTAGCAGTGACTATGTCTGGATCTATCTGCTT 9650
PVYN|X97895 ATCGATTAGGTGATGTTGCGA-TTGTCTAGCAGTGACTATGTCTGGATTTAGTTACTT 9648
***** ** * * * * * * * * * * * * * * * * * * * *

PVYO|AJ585195 GGGTGGTGTGTTGATTTTCGTCATAGCAGTGACTG----- 9651
PVYO|AJ585196 GGGTGGTGTGTTGATTTTCGTCATAACAGTGACTGTAACTTCAATCAGGAGAC 9698
PVYO|AJ890349 GGGTGGTGTGTTGATTTTCGTCATAACAGTGACTGTAACTTCAATCAGGAGAC 9699
PVYO|DQ157178 GGGTGGTGTGTTGATTTTCGTCATAACAGTGACTGTAACTTCAATCAGGAGAC 9670
PVYO|EF558545 GGGTGGTGTGTTGATTTTCGTCATAACAGTGACTGTAACTTCAATCAGGAGAC 9698
PVYO|HQ912864 GGGTGGTGTGTTGATTTTCGTCATAACAGTGACTGTAACTTCAATCAGGAGAC 9698
PVYO|HQ912865 GGGTGGTGTGTTGATTTTCGTCATAACAGTGACTGTAACTTCAATCAGGAGAC 9698
PVYO|HQ912888 GGGTGGTGTGTTGATTTTCGTCATAACAGTGACTGTAACTTCAATCAGGAGAC 9667
PVYO|U09509 GGGTGTGTTGTTGATTTTCGTCATAACAGTGACTGTAACTTCAATCAGGAGAC 9698
PVYO|X12456 GGGTGTGCTGTGATTCTGTCTATAGCAGTGACTGTAACTTCAATCAGGAGAC 9704
PVYN|AJ585197 GGGTGTGCTGTGATTCTGTCTATAGCAGTGACTG----- 9653
PVYN|AJ585342 GGGTGGTGTGTTGATTTTCGTCATAACAGTGACTG----- 9647
PVYN|AJ889866 GGGTGGTGTGTTGATTTTCGTCATAACAGTGACTGTAACTTCAATCAGGAG-- 9697
PVYN|AY166866 GGGTGTGCTGTGATTCTGTCTATAGCAGTGGCTGTAACTTCAATCAGGAGAC 9700
PVYN|AY884982 GGGTGGTGTGTTGATTTTCGTCATAACAGTGACTGTAACTTCAATCAGGAGAC 9702
PVYN|AY884983 GGGTGTGTTGTTGATTTCTGTCTATAGCAGTGACTGTAACTTCAATCAGGAGAC 9700
PVYN|FJ666337 GGGTGGTGTGTTGATTTTCGTCATA----- 9583
PVYN|HQ912869 GGGTGGTGTGTTGATTTTCGTCATAACAGTGACTGTAACTTCAATCAGGAGAC 9671
PVYN|M95491 GGGTGGTGTGTTGATTTTCGTCATAACAGTGACTGTAACTTCAATCAGGAGAC 9703
PVYN|X97895 GGGTGTGCTGTGATTCTGTCTATAGCAGTGACTGTAACTTCAATCAGGAGAC 9701
***** ** * * * * *

```

Example 2.2. Multiple sequence alignment of PVY strains that belong to different serotypes.

The processing of such large nucleotide sequences like viral genomes is required long computation time, therefore, the Clustal Omega user can to obtain the alignment by email after few minutes or hours. The computation time usually depends from number of strains and their sequence length. The user should to save the alignment as text file with *.txt extension.

Step 5. Masking serotype-specific mutations.

Here, we again suggest to use the MorphoCatcher service (<http://morphocatcher.ru>) for the MSA processing. To visualize genome region with high density of serotype-specific mutations, the sliding-window plot of average mutation index should be obtained (**Figure 2.2**).




Figure 2.2. The alignment-based average mutation index plot for two PVY serotypes.

Step 6. Designing specific primers.

Finally, the user should to upload corresponding text file from the MorphoCatcher output (see files from the "Input for PrimerExplorer" column) as the input file to the PrimerExplorer service (<http://primerexplorer.jp/e>) for designing of serotype-specific primer set.

Here, the fragment of primer set list is shown (**Figure 2.3**) to demonstrate how different serotype-specific mutations can be covered by the primer ends of various primer sets.




Figure 2.3. The PrimerExplorer output page with full view of designed primer sets. Specific mutations for PVY strains of O-serotype are highlighted by red color. The same mutations can be covered by oligos from different primer sets.

Credits

We are grateful to the "Start Bootstrap" service (<https://startbootstrap.com>) for providing a free "Agency" template of HTML5 page. Some graphics for web design was provided by contributors of the "Noun Project" collection of icons (<https://thenounproject.com>) — we thank Danil Polshin ("Abstract" icon), Wes Breazell ("DNA" icon), Alexander Blagochevsky ("Pedigree Chart" icon), Ervin Bolat ("Text Editor" icon), Cris Dobbins ("Sliders" icon), and ProSymbols ("Idea" icon).

References

1. PubMed ID: 25684775;
2. PubMed ID: 29119360;
3. PubMed ID: 21675922.